

Michelson's Interferometer


An interferometer is an instrument in which the phenomenon of interference is used to precise measurements of wavelengths or distances.

Principle

- 1) A beam of light from extended source splits into two coherent beams of equal intensities by partial reflection or refraction.
- 2) These beams travel in two mutually perpendicular paths & meet after reflection from plane mirrors.
- 3) The beams overlap on each other to form interference fringes.
- 4) If mirrors have circular small angle inclination \rightarrow straight

Construction

- 1) It has two plane mirrors M_1 & M_2 .
- 2) Beam splitter $G_1 \rightarrow$ partially silvered plane II glass plate
 \rightarrow transmits half the incident light & reflects back the rest.
- 3) Compensating plate $G_2 \rightarrow$ simple plane II glass plate having thickness same as G_1 .



- 4) G_1, G_2 are \perp to each other & inclined at an angle 45° to mirror M_2 .
 - 5) M_1 is mounted on the carriage & can be moved ~~to~~ exactly 10 to itself with micrometer screws.
 - 6) M_1 can be moved with help of graduated drum attached to the screw.
 - 7) M_1, G_1, M_2 can be made ~~to~~ exactly 1 mm with help of fine screws attached to it. & interference bands are observed in the fields of view of telescope T .
- $AB = 1$

Path difference calc.

Path of wave 1

$$AM_1 AT + \frac{\lambda}{2} + \frac{\lambda}{2} = AM_1 AT + \lambda$$

$$= AM_1 A + \lambda = 2x_1 + \lambda \rightarrow \textcircled{1}$$

Path of wave 2

$$AM_2 AT + \frac{\lambda}{2} + \frac{\lambda}{2} = AM_2 AT + \lambda$$

$$= AM_2 A + \lambda = 2x_2 + \lambda \rightarrow \textcircled{2}$$

\therefore path diff b/w 2 waves eq^n

$$S = 2(x_2 - x_1)$$

when it is viewed like a thin film, then the reflected waves are

$$AB = L$$

causing a path diff
since $n_2 - n_1 = d$

Thickness of the thin film i.e. distance b/w mirror M_1 & virtual M_2'

$$\therefore S = 2d$$

When the film is viewed at an angle of θ the path diff can be written as

$$S = 2d \cos \theta + \frac{\lambda}{2}$$

case(i) $d=0$

$$\therefore S = \frac{\lambda}{2}$$

for constructive $\Rightarrow S = n\lambda$

for destructive $\Rightarrow S = (2n+1)\frac{\lambda}{2}$

case(ii) distance b/w $M_1, M_2' = \frac{1}{u}$

$$S = \frac{\lambda}{2} + \left(\frac{1}{u} + \frac{1}{u} \right) = \frac{\lambda}{2} + \frac{\lambda}{2} = \lambda$$

\therefore This should result in bright fringe

case(iii) increase distance with $\frac{1}{n_L}$

$$S = \frac{d}{2} + \left(\frac{d}{n} + \frac{d}{n} \right) + \left(\frac{d}{n} + \frac{d}{n} \right)$$

$$S = \frac{d}{2} + \left(\frac{d}{2} \right) + \left(\frac{d}{2} \right)$$

$$S = \frac{3d}{2}$$

which results in dark fringe.

Applications of Interferometer

- 1] To determine wavelength of incident light.
- A monochromatic light is placed at the position of the source 'S' whose wavelength needs to be determined. Let d_1 be the distance b/w the source and mirror M_1 , whose reading distance of micrometer of position M_1 . Now change the position of M_1 to get d_2 . The no. of fringes disappearing during the same time is taken as N . Therefore wavelength of the incident light can be determined by

$$\lambda = \frac{2(d_2 - d_1)}{N}$$

- 2) Difference b/w the wavelengths of two spectral lines which are close to each other.

Determining the avg wavelength of the two spectral lines, the diff. b/w wavelengths is given of the two lines is given by

$$\Delta\lambda = \frac{\lambda^2}{2d}$$

d = distance b/w M_1, C, M_2'

- 3) Thickness of a thin transparent sheet.

For 'N' no. of fringes displaced in a fringe system, when insertions of thin transparent sheets takes place in the path of one of the interfering beams then thickness is given by

$$t = \frac{N\lambda}{2(\mu-1)}$$

μ = RI of the sheet.